

Location \_\_\_\_\_

### IMLS Survey - Textiles & Costumes

Acc# \_\_\_\_\_  Textile  Costume Description: \_\_\_\_\_

---

**Condition:**       Stable                       Needs rehousing                       Needs first aid

30 min - rolled                       30 min - padded hanger - 20th c.  
 30 min - boxed                       60 min - padded hanger - 18th/19th c.  
 90 min - custom mount                       30 min - padded hanger - childs  
 60 min - custom mount                       30 min - padded hanger - skirt/pant  
 30 min - custom mount  
 60 min - custom box

---

**Proposed Storage:**

**Drawer**       Flat              L: \_\_\_\_\_      W: \_\_\_\_\_      H: \_\_\_\_\_  
 **Tray**  
 **Either**       Folded              L: \_\_\_\_\_      W: \_\_\_\_\_      H: \_\_\_\_\_  
 Rolled              L: \_\_\_\_\_      D: \_\_\_\_\_      Tube size:  1 1/2"  2"  3"  4"

**Box**               30" x 18" x 6"     30" x 22" x 4"     30" x 24" x 6"     40" x 18" x 6"     60" x 18" x 6"  
 60" x 36" x 5"     shoe 12" x 8" x 5"     boot 14" x 7" x 12 1/2"

whole     1/2     1/3     1/4     1/6     1/8

**Custom Box**    L: \_\_\_\_\_      W: \_\_\_\_\_      H: \_\_\_\_\_

**Shelf (fixed)**    L: \_\_\_\_\_      W: \_\_\_\_\_      H: \_\_\_\_\_

custom mount

**Rolled**              L: \_\_\_\_\_      Diam.: \_\_\_\_\_      Tube size:  1 1/2"  2"  3"  4"

**Rack hanging**    H: \_\_\_\_\_      W: \_\_\_\_\_      D: \_\_\_\_\_  
 S-hooks               Velcro

**Needs special handling during move**